

Part 1 – Preparation of Analyte Sample**Partner #1**

1. Obtain the unknown steel sample-bottle corresponding to the number assigned to you and your partner by your instructor. *This is not the same as your locker number.*
2. Using a metal scoopula, accurately weigh into a tared weigh-boat **0.20 to 0.30 g** of the powdered steel at an analytical balance. Seal and return the sample-bottle.
3. Carefully add a few drops of distilled water to the steel powder to moisten it.¹

Quantitatively transfer the steel into a clean, but not necessarily dry, 250 mL Erlenmeyer flask, using the smallest amount of distilled water as possible, as described below:

Begin by squeezing-together two sides of the weigh-boat to create a lip at one end. Carefully tilt the lip of the weigh-boat over the mouth of the Erlenmeyer flask, allowing the wet powder to roll into the flask. Then use a gentle spray of a very small amount of distilled water to rinse any remaining steel from the weigh boat into the flask. It is normal for a small brown stain to remain in the weigh boat; this is not steel.

4. **WEARING PROTECTIVE GLOVES**, add **40 mL** of **5 M nitric acid (HNO₃)** directly to the wet steel powder in the 250 mL Erlenmeyer flask. *A pump-bottle of the acid will be used.*
5. **AT A FUMEHOOD** place the Erlenmeyer flask on a hotplate (already set up for you by the lab technician). Turn the heat-setting on to almost maximum and wait for the liquid in the flask to reach a boil. Once boiling, **continue to heat for 10 minutes**, noting that:
 - after about 3 minutes of boiling, ensure that no metal-particles are stuck on the upper-walls of the flask. If some are, rinse them down using a small amount of water.
 - **the sample must not be allowed to boil dry as the flask will shatter.** If necessary, make up any lost volume by periodic addition of small quantities of distilled water.
 - a toxic brown gas of NO₂ is produced by this reaction. **The fumehood's window-sash must be lowered to the point labeled by a red-arrow on the side of the fumehood's window-frame to ensure no inhalation of this gas.**
 - all of the steel in the sample will have dissolved after 10 minutes, resulting in a clear yellow or orange solution due to the presence of Fe³⁺. Manganese is present as the colourless cation Mn²⁺. Some black particles of carbon may remain undissolved.
6. Sometime during the prior 10 minute boiling-phase, weigh out approximately **1 g** of ammonium persulphate, **(NH₄)₂S₂O₈**, into a plastic weigh-boat using a top-loading balance. Do **not** add it to your flask yet. Set this white powder aside at the fumehood, ready for use later on in step #8. *The mass weighed is not recorded on your report sheet.*

¹ The electrostatic steel particles tend to fly away when transferred dry, hence the need to moisten them.

- When the steel sample has boiled for 10 minutes, use your crucible-tongs to move the hot flask onto the silver heat-resistant mat that sits next to the hotplate in the fumehood. *Never place hot glassware on the plastic liner in the fumehood, as it will instantly melt!*
- While still at the fumehood, use a scoopula to **add approximately half** of the previously-weighed $(\text{NH}_4)_2\text{S}_2\text{O}_8$ into the hot flask. Using crucible-tongs, carefully give the hot flask a quick swirl to mix the contents. *The ensuing chemical reaction causes the solution to bubble or froth-up as the excess nitric acid and nitric oxides are destroyed.*

Once the chemical reaction has calmed, add the remaining **second-half** of ammonium persulphate, $(\text{NH}_4)_2\text{S}_2\text{O}_8$, to the flask and swirl it as before.

Return the flask to the hot-plate and boil gently for 10 minutes. *Check your flask periodically to ensure it doesn't boil dry; add distilled water if needed to prevent this.*

- Once the boiling of step #8 is over, use crucible-tongs to set the hot flask onto the heat-resistant mat. **Dilute the sample to ~ 100 mL with distilled water**, using the lines on the side of the flask for guidance.
- Let the flask cool down a bit until you can hold the neck safely with your gloved-hands.

Take the flask to the fumehood where a pump-bottle of **concentrated (85%) phosphoric acid**, H_3PO_4 is located. Pump **10 mL** of the H_3PO_4 directly into your flask. *The acid forms a complex with Fe(III), thereby allowing the purple-colour of Mn^{7+} in step #12 to be visible.*

- Weigh out approximately **0.5 g** of potassium periodate, KIO_4 , into a plastic weigh-boat at a top-loading balance. *This chemical is also a white powder; the mass is not recorded.*
- At the fumehood, pour all of the KIO_4 powder into the flask, swirl to mix, and then return the flask to the hot-plate. Wait for the contents to reach a boil. Once at this temperature, **boil the sample for just 1 minute.**

Reduce the heat-setting on the hotplate to the minimum setting and leave the flask there for 10 minutes to heat gently. *The solution will very gradually change to a purple colour as the manganese oxidizes from Mn^{2+} to Mn^{7+} , in the form of MnO_4^- .*

- At a sink, **cool the purple solution to room temperature** by running a stream of cold water against the outside of the Erlenmeyer flask, with swirling of the contents.
- Refer to Figure 2.3** in the prior Titration Experiment # 2.2 (Part 1 – KHP). Apply all of those techniques to transfer the cooled purple MnO_4^- analyte solution quantitatively from the Erlenmeyer flask into a rinsed 250.00 mL volumetric flask.

Dilute carefully to the 250 mL calibration line with distilled water. **DO NOT OVERSHOOT!** Cap the flask and mix well to ensure a uniform solution. The MnO_4^- analyte is now ready for colorimetric analysis in Part 3.