


for students at locker # 14



Trial 1: Hydrolysis of tBuCl at Room Temperature [tBuCl] vs time; (t = 0 - 300 s)



MN Station # 2 GRAPHS p. 3



These are the First and Second-Order plots for Trial #1.

Note that data is not being evaluated over the entire 300 seconds because as we approach t = 300 seconds, [tBuCl] approaches 0.0 M. The first-order operation 'ln(0)' is undefined, while the second-order operation '1 / 0' yields infinity. Plotting such results would yield useless graphs.


These graphs might look slightly different to what you are accustomed to seeing in a textbook, as those graphs usually only show 5 or 6 data-points, and then a line-of-best-fit drawn through them. With our MeasureNet monitored experiments, we have gathered hundreds of data-points, as a pH-reading was taken every 0.5 seconds. The solid line on each graph arises from all of these data-points.


The order of the reaction with respect to tBuCl can be determined by whichever of these two graphs yields a straight solid-line.

For both graphs, the dashed-line is the Excel-generated 'line-of-best-fit'. The equations shown are of the form 'y = mx + b', and the R^2 value is a measure of how well the data fits the line, where an R^2 of 1.0 indicates a perfect fit.

From the integrated rate-laws, the rate-constant, k, can be calculated from the slope (m) of the graphs:

First-order: $k = - slope \quad units \ of \ s^{-1}$ Second-order: $k = slope \quad units \ of \ M^{-1} \ s^{-1}$ MN Station # 2 GRAPHS p. 4





This graph shows the first-order plot for Trials #1 - #3, but only over the time-domain of 0 - 2 seconds so that the very initial rate can be calculated, when the reaction is at it's fastest.


The line-of-best-fit is shown for all three trials, and the table shows the slope and R² values for each line.

	y = mx + b	
trial #1:	$y = -0.175 \times -5.17$	$R^2 = 0.9978$
trial #2:	$y = -0.498 \times -5.13$	$R^2 = 0.9798$
trial #3:	y = -0.0838 x - 5.17	$R^2 = 0.9931$

As discussed on Graph #3, the rate-constant, k, is the negative of the slope of the best-fit line.

- Convert the slope-values to initial rate-constants, k, for all three trials and record these on your Lab Report Sheet (in units of s^{-1}).
- Use the appopriate equation from the lab manual to calculate the theoretical first half-lifes.

You will now need to use the MS-Excel workbook called **KIN2_Arrhenius_Plotter** to generate an Arrhenius Plot. Your instructor will have already sent this Excel file to your Capilano University email account. Once you have obtained the file, open it using MS-Excel and follow the on-screen instructions.



This graph is similar to Graph #2, except that here, the y-axis is a measure of what percent of unreacted tBuCl remains in the beaker. The conversion from a [tBuCl] y-axis to a % tBuCl remaining y-axis is simple. Each data-point (in units of tBuCl M) was simply divided by the initial concentration of tBuCl at t = 0 seconds, and then multiplied by 100. Note that only the first 60 seconds of data have been plotted, so as to magnify the first portion of the x-axis.

Perform on this page the Graphical Determination of Half-Lives:

- find the point on each curve where it intersects with the 50% tBuCl Remaining line.
- drop a vertical-line down from each point to the x-axis.
- estimate, to the nearest 1 s, the time where each dropped-line intersects the x-axis and then record these on your Report Sheet as the First Half-Life.

Be sure to show all of this work directly on this graph page.