

RAW DATA: Ideal Gas Law**Section L2 - TWO**

Please record the necessary data for your locker number.

Locker	volume flask (mL)	initial mass (g)	final mass (g)	temp (°C)	pressure (kPa)
A202	145.0	75.3989	75.6678	84.2	102.0
A204	157.0	92.7017	93.1284	90.0	102.4
A207	158.9	93.0668	93.2525	85.8	102.7
A211	156.0	92.2334	92.5206	88.5	101.4
A213	150.0	95.4375	95.6203	79.1	102.2
A217	148.0	93.6871	94.0091	80.7	101.1
A219	148.5	80.8740	81.3457	88.9	101.8
A220	151.5	93.2042	93.3536	90.0	102.1
B202	155.0	93.1126	93.3916	88.5	102.4
B204	147.1	96.6906	97.0862	84.2	101.9
B207	148.0	76.0903	76.2589	88.0	101.7
B211	149.0	96.4890	96.7966	87.5	101.5
B213	149.5	79.9618	80.4261	84.3	102.0
B217	157.5	93.6217	93.7846	89.0	102.3
B219	153.3	94.7979	95.0786	93.2	101.1
B220	151.9	76.5138	76.9290	92.5	101.5