

RAW DATA: Ideal Gas Law**Section L3 - THREE**

Please record the necessary data for your locker number.

Locker	volume flask (mL)	initial mass (g)	final mass (g)	temp (°C)	pressure (kPa)
A302	146.0	80.1383	80.3070	91.5	102.2
A304	147.5	82.2834	82.5942	89.0	102.1
A307	148.0	97.1710	97.6339	87.5	101.7
A311	150.0	94.8198	94.9703	89.8	101.6
A313	145.0	75.3989	75.6678	84.2	102.0
A317	157.0	92.7017	93.1284	90.0	102.4
A319	158.9	93.0668	93.2525	85.8	102.7
A320	156.0	92.2334	92.5206	88.5	101.4
B302	150.0	95.4375	95.6203	79.1	102.2
B304	148.0	93.6871	94.0091	80.7	101.1
B307	148.5	80.8740	81.3457	88.9	101.8
B311	151.5	93.2042	93.3536	90.0	102.1
B313	155.0	93.1126	93.3916	88.5	102.4
B317	147.1	96.6906	97.0862	84.2	101.9
B319	148.0	76.0903	76.2589	88.0	101.7
B320	149.0	96.4890	96.7966	87.5	101.5