

RAW DATA: Ideal Gas Law**Section L5 - FIVE**

Please record the necessary data for your locker number.

Locker	volume flask (mL)	initial mass (g)	final mass (g)	temp (°C)	pressure (kPa)
A502	155.0	93.1126	93.3916	88.5	102.4
A504	147.1	96.6906	97.0862	84.2	101.9
A507	148.0	76.0903	76.2589	88.0	101.7
A511	149.0	96.4890	96.7966	87.5	101.5
A513	149.5	79.9618	80.4261	84.3	102.0
A517	157.5	93.6217	93.7846	89.0	102.3
A519	153.3	94.7979	95.0786	93.2	101.1
A520	151.9	76.5138	76.9290	92.5	101.5
B502	146.0	80.1383	80.3070	91.5	102.2
B504	147.5	82.2834	82.5942	89.0	102.1
B507	148.0	97.1710	97.6339	87.5	101.7
B511	150.0	94.8198	94.9703	89.8	101.6
B513	145.0	75.3989	75.6678	84.2	102.0
B517	157.0	92.7017	93.1284	90.0	102.4
B519	158.9	93.0668	93.2525	85.8	102.7
B520	156.0	92.2334	92.5206	88.5	101.4